

GD/F30S Gynecological Training Models

This Gynecological Training Model is specifically designed for professional training.

The various functions of these model are as follows :-

1. Palpation of normal and pregnant uterus
 2. Bi-manual pelvic and tri-manual pelvic examination
 3. Vaginal examination
 4. Visual recognition of normal and abnormal cervixes
 5. Contraceptive sponge insertion and removal
- Observation of uterus, ovary, fallopian tubes and ligament teres.




GD/FW2 Maternity Examination Model

This Advanced Maternity Examination Model is specifically designed for professional training.

The various functions of these models are as follows :-

1. Palpation, embryo heart sound auscultation, pelvis measurement and mammary nursing training
2. Inflating uterus through ball to make model lifelike
3. Adjustable heart sound speed, volume
4. Heart sound can be heard not only from stethoscope but also from microphone in front board of the model



GD/F3B Vulva Suturing Training Simulator

The models demonstrate the wound in different positions in pudendum. It's made of soft materials, can be used many times in suturing training.

