


## GD/H100S Nursing Training Manikin (Unisex)

This full-body lifelike manikin is specifically designed for professional nursing training. This high quality manikin will provide wide range of exercises possible including both first aid and emergency situations.



### Features:

- |   |  |
|---|--|
| <ol style="list-style-type: none"> <li>1. Hair and face washing</li> <li>2. Eye and ear washing and administering</li> <li>3. Mouth cavity and artificial teeth care</li> <li>4. Endotracheal intubation</li> <li>5. Sputum suction</li> <li>6. Oxygen inhaling</li> <li>7. Oral and nasal feeding</li> <li>8. Gastrolavage</li> <li>9. Main organs in thoracic cavity</li> <li>10. Venipuncture, injection, blood transfusion (Arm)</li> <li>11. Deltoid subcutaneous injection</li> <li>12. Vastus lateralis injection</li> </ol> | <ol style="list-style-type: none"> <li>13. Thoracic cavity, abdominal cavity, liver, bone marrow and lumbar puncture</li> <li>14. Enema</li> <li>15. Female urethral Catheterization</li> <li>16. Male urethral catheterization</li> <li>17. Female bladder irrigation</li> <li>18. Male bladder irrigation</li> <li>19. Ostomy</li> <li>20. Buttocks intramuscular injection</li> <li>21. Main organs in abdominal cavity</li> <li>22. Holistic nursing care: sponge bath, replacing clothes</li> </ol> |
|---|--|

## GD/H2 Multi-Functional Basic Nursing Manikin (Female)

This Nursing Manikin can demonstrate the following Basic Nursing Skills.



### Features:

- |   |  |
|---|--|
| <ol style="list-style-type: none"> <li>1. Hair and face washing, hair combing and sponge bath</li> <li>2. Mouth cavity care</li> <li>3. Endotracheal intubation</li> <li>4. Oxygen inhaling</li> <li>5. Nasal feeding</li> <li>6. Gastrolavage</li> <li>7. Intracardiac injection</li> <li>8. Closed chest cardiac resuscitation</li> <li>9. Pneumothorax</li> <li>10. Breast Care</li> <li>11. Liver puncture</li> </ol> | <ol style="list-style-type: none"> <li>12. Chest drainage</li> <li>13. Abdominal drainage</li> <li>14. Bone marrow puncture</li> <li>15. Lumbar puncture</li> <li>16. Deltoid subcutaneous injection</li> <li>17. Venipuncture injection</li> <li>18. Blood transfusion (Arm)</li> <li>19. Vastus lateralis injection</li> <li>20. Female urethral catheterization</li> <li>21. Female bladder irrigation</li> <li>22. Buttocks intramuscular injection</li> </ol> |
|---|--|