

GD/LV2 Advanced Surgical Suture Leg

This model of Advanced Surgical Suturing Leg allows the students to practice cut, suturing, stitches release, enswathing, etc. The elasticity and flexibility of skin can make it for practice hundreds of times. There are some initial wounds on the arm and if all the initial wounds have been sutured, people can cut the arm in other place for new training of suturing.



GD/J5S Airway Intubation Model

Features:

1. Accurate anatomical marks
2. Oral or nasal cavity intubation
3. Wrong operation, stomach inflation or laryngoscope compressing the teeth will activate warning prompting
4. One proper pupil and another mydriasis



GD/H10C Adult Intramuscular Injection Training Mannequin



Features:

1. Structure of left gluteal: osteal marks, gluteal muscles, ischiadic nerve and veins
2. Left gluteal can be taken off. So it is easy to observe internal structure and make sure position of ischiadic nerve and veins.
3. Gluteal intramuscular injection training

NC-6028 Buttock Inter Muscular Injection Training Model



A well-Shaped buttocks of an adult, made of special material, correct anatomical landmarks.

- Adult buttock, precise anatomical structure
- Both sides are suitable for injection practice, realistic texture
- Electronic monitoring, indicating injection position, depth by LED lights and sound alert.
- Gluteus maximus injection, allow injecting "portion" convenient drain outlet
- Gluteus medius, gluteus minimus